

BIVONA® silicone tracheostomy tube range

The softer approach to tracheostomy


















AIRWAY MANAGEMENT

For over 70 years, the name of Portex® has been synonymous with the values of quality, reliability and innovation. These core values are as strong as they have ever been, with constant investment in quality systems and research and development.

Our range of Bivona® tubes embrace these core values to provide a range of high quality silicone products that ultimately give your patients a softer approach to tracheostomy.


Contents

Page

Bivona® Silicone Tubes		4-5
Bivona®: For Choice, Quality & Design		6-7
Bivona® TTS™		8-9
Bivona® Aire Cuf®		10-11
Bivona® Fome-Cuf®		12-13
Bivona® Cuffless		14-15
Bivona® Adjustable Hyperflex™		16-17
Bivona® Fixed Flange Hyperflex™		18-19
Bivona® Flextend™		20-21
Bivona® FlexTend™ TTS™		22-23
Bivona® Sleep Apnea Tubes		24-25
Bivona® Laryngectomy Tubes		24-25
Bivona® Epistaxis Catheter		26
Bivona® Customised		27
Accessories		28

Silicone is a material with many ideal qualities, offering important advantages over conventional tube materials:

- Softness and flexibility
 - For comfort and ease of use
- Biocompatibility
 - Hydrophobic and hypoallergenic, tissue-friendly and easy to clean
- Durability
 - Robust for reliable reprocessing, withstands repeated reprocessing up to ten times for an adult tube and five times for a neonatal or paediatric tube.



Bivona[®] silicone tubes add advanced design features to the benefits of silicone. Available in a range of tube configuration options to suit each age, size and application requirement. Customised options are also available.

For quality.

Quality is our passion in everything we do. You can be assured of the highest possible standards when using our range of Bivona® tracheostomy tubes.

For choice.

Combining the softness and flexibility of silicone, with an extensive range of size and shape options for all ages. We are unique in offering a silicone tube that has the additional benefits of an Inner Cannula, making the Bivona® range a perfect solution for patients whether at home or in a hospital.

For design.


Integrating unique features that may help to aid patient comfort and ease of use.

For useability and safety.

Hydrophobic material for easy cleaning, free from DEHP, latex and plasticizers.

For cost-efficiency.

Bivona® tubes can be cleaned and reprocessed up to ten times for adult sizes and five times for neonatal and paediatric sizes.



Bivona®

TTS™

Tight to Shaft

The TTS™ cuff design offers the benefits of a cuff with the profile of a cuffless tube. When the silicone cuff is totally deflated it adds no significant dimension to the outer diameter of the tube's shaft, yet it can be partially or completely inflated for varied airway management needs.

1



670180

2



67N025

3



67P040

4





67SN025

5



67SP040



= Inner Cannula Compatible



= Recommended for hospital use



= Recommended for home use



1

Bivona® TTS™ Cuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
670150	5.0 mm	5.0 mm	7.4 mm	60.0 mm	1
670160	6.0 mm	6.0 mm	8.8 mm	70.0 mm	1
670170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
670175	7.5 mm	7.5 mm	10.5 mm	80.0 mm	1
670180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
670185	8.5 mm	8.5 mm	11.8 mm	88.0 mm	1
670190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1



2

Bivona® TTS™ Cuffed Neonatal “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
67N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
67N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
67N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

3

Bivona® TTS™ Cuffed Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
67P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
67P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
67P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
67P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
67P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
67P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

4

Bivona® TTS™ Cuffed Neonatal “Straight” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67SN025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
67SN030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
67SN035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
67SN040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

5

Bivona® TTS™ Cuffed Paediatric “Straight” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67SP025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
67SP030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
67SP035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
67SP040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
67SP045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
67SP050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
67SP055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

Bivona[®]

Aire-Cuf[®]

A flexible, adjustable tube with an air-filled cuff. Neonatal and paediatric are reinforced to give a high level of kink and crush resistance.

1



750160

2



755180

3



65N025

4



65P055

5




65SN035

6



65SP050



= Inner Cannula Compatible



= Recommended for hospital use



= Recommended for home use



1

Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
750150	5.0 mm	5.0 mm	7.4 mm	60.0 mm	1
750160	6.0 mm	6.0 mm	8.8 mm	70.0 mm	1
750170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
750180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
750190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
750195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



2

Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes with Talk Attachment

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
755150	5.0 mm	5.0 mm	7.4 mm	60.0 mm	1
755160	6.0 mm	6.0 mm	8.8 mm	70.0 mm	1
755170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
755180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
755190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
755195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1

3

Bivona® Aire-Cuf® Neonatal “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
65N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
65N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
65N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

4

Bivona® Aire-Cuf® Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
65P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
65P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
65P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
65P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
65P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
65P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

5

Bivona® Aire-Cuf® Neonatal “Straight” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65SN025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
65SN030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
65SN035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
65SN040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

6

Bivona® Aire-Cuf® Paediatric “Straight” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65SP025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
65SP030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
65SP035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
65SP040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
65SP045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
65SP050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
65SP055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

Bivona[®]

Fome-Cuf[®]

This innovative cuff expansion system expands to provide a tracheal seal to prevent leaks and protect against aspiration. The Fome-Cuf[®] reduces the risk of tracheal dilation by automatically adjusting to changes in the shape of the patient's trachea.

1



850180

2



855180

3



85N030

4



85P035



= Recommended for hospital use



= Recommended for home use



1

Bivona® Fome-Cuf® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
850150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
850160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
850170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
850180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
850190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
850195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1

2

Bivona® Fome-Cuf® Tracheostomy Tubes with Talk Attachment

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
855150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
855160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
855170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
855180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
855190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1

3

Bivona® Fome-Cuf® Neonatal "V" Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
85N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
85N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
85N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

4

Bivona® Fome-Cuf® Paediatric "V" Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
85P025	2.5mm	2.5mm	4.0 mm	38.0 mm	1
85P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
85P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
85P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
85P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
85P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
85P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

Bivona®

Cuffless

Bivona® uncuffed tracheostomy tubes with their all-silicone construction are intended to provide maximum comfort. The silicone tube material is hydrophobic and hypoallergenic, allowing for a tube that is tissue-friendly and easy to clean.

1



60A160

2



60N040

3



60P040

4





60SN025

5



60SP045



= Inner Cannula Compatible



= Recommended for hospital use



= Recommended for home use



1

Bivona® Uncuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60A150	5.0 mm	5.0 mm	7.4 mm	60.0 mm	1
60A160	6.0 mm	6.0 mm	8.8 mm	70.0 mm	1
60A170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
60A180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
60A190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1



2

Bivona® Uncuffed Neonatal "V" Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60N025	2.5mm	2.5mm	4.0 mm	30.0 mm	1
60N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
60N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
60N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

3

Bivona® Uncuffed Paediatric "V" Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
60P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
60P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
60P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
60P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
60P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
60P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

4

Bivona® Uncuffed Neonatal Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60SN025	2.5 mm	2.5mm	4.0 mm	30.0 mm	1
60SN030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
60SN035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
60SN040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1

5

Bivona® Uncuffed Paediatric Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60SP025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
60SP030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
60SP035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
60SP040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
60SP045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
60SP050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
60SP055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1

Bivona[®]

Adjustable

Hyperflex[™]

A one-size-fits-all tube that can be instantly customised to fit your patient. It can be adjusted for horizontal and vertical shaft drop to accommodate unusual anatomy or pathology. Flexible wire-reinforced silicone shafts provide kink and crush resistance.

1



60HA35

2




67HA70

3



75HA70



= Recommended for hospital use



= Recommended for home use



1

Bivona® Uncuffed Adjustable Neck Flange Hyperflex™ Paediatric“V” Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60HA25	2.5 mm	2.5 mm	4.0 mm	56.0 mm	1
60HA30	3.0 mm	3.0 mm	4.7 mm	61.0 mm	1
60HA35	3.5 mm	3.5 mm	5.3 mm	66.0 mm	1
60HA40	4.0 mm	4.0 mm	6.0 mm	71.0 mm	1
60HA45	4.5 mm	4.5 mm	6.7 mm	76.0 mm	1
60HA50	5.0 mm	5.0 mm	7.3 mm	81.0 mm	1
60HA55	5.5 mm	5.5 mm	8.0 mm	86.0 mm	1

2

Bivona® TTS™ Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67HA60	6.0 mm	6.0 mm	9.2 mm	110.0 mm	1
67HA70	7.0 mm	7.0 mm	10.6 mm	120.0 mm	1
67HA80	8.0 mm	8.0 mm	11.7 mm	130.0 mm	1
67HA90	9.0 mm	9.0 mm	12.9 mm	140.0 mm	1

3

Bivona® Aire-Cuf® Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
75HA60	6.0 mm	6.0 mm	9.2 mm	110.0 mm	1
75HA70	7.0 mm	7.0 mm	10.6 mm	120.0 mm	1
75HA80	8.0 mm	8.0 mm	11.7 mm	130.0 mm	1
75HA90	9.0 mm	9.0 mm	12.9 mm	140.0 mm	1

Bivona[®]

Fixed Flange

Hyperflex[™]

With a longer than standard shaft, this tube provides both extra vertical and horizontal length for a more permanent solution. Flexible wire-reinforced silicone shaft provides improved kink and crush resistance.

1



60AFHXL70

2




67FHXL70

3



75FHXL60



= Recommended for hospital use



= Recommended for home use



1

Bivona® Uncuffed Extra Length Fixed Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60AFHXL60	6.0 mm	6.0 mm	9.2 mm	110.0 mm	1
60AFHXL70	7.0 mm	7.0 mm	10.6 mm	120.0 mm	1
60AFHXL80	8.0 mm	8.0 mm	11.7 mm	130.0 mm	1
60AFHXL90	9.0 mm	9.0 mm	12.9 mm	140.0 mm	1

2

Bivona® TTS™ Extra Length Fixed Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67FHXL60	6.0 mm	6.0 mm	9.2 mm	100.0 mm	1
67FHXL70	7.0 mm	7.0 mm	10.6 mm	110.0 mm	1
67FHXL80	8.0 mm	8.0 mm	11.7 mm	120.0 mm	1
67FHXL90	9.0 mm	9.0 mm	12.9 mm	130.0 mm	1

3

Bivona® Cuffed Extra Length Fixed Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
75FHXL60	6.0 mm	6.0 mm	9.2 mm	110.0 mm	1
75FHXL70	7.0 mm	7.0 mm	10.6 mm	120.0 mm	1
75FHXL80	8.0 mm	8.0 mm	11.7 mm	130.0 mm	1
75FHXL90	9.0 mm	9.0 mm	12.9 mm	140.0 mm	1

Bivona® FlexTend™

Reinforced proximal and distal sections provide a high level of kink and crush resistance. Extended tubes may help minimise chin abrasions on “no-neck” patients.

1



60NFP30

2



60NFP35

3



60PFS55

4



60PFSS045

5




60PFP50

6



60PFP50



= Recommended
for hospital use



= Recommended
for home use



1

Bivona® FlexTend™ Uncuffed Neonatal “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60NFP25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	30.0 mm	1
60NFP30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	32.0 mm	1
60NFP35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	34.0 mm	1
60NFP40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	36.0 mm	1

2

Bivona® FlexTend™ Uncuffed Neonatal Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60NFPS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	30.0 mm	1
60NFPS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	32.0 mm	1
60NFPS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	34.0 mm	1
60NFPS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	36.0 mm	1

3

Bivona® FlexTend™ Uncuffed Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	38.0 mm	1
60PFS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	39.0 mm	1
60PFS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	40.0 mm	1
60PFS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	41.0 mm	1
60PFS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	42.0 mm	1
60PFS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	44.0 mm	1
60PFS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	46.0 mm	1

4

Bivona® FlexTend™ Uncuffed Paediatric Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFSS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	38.0 mm	1
60PFSS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	39.0 mm	1
60PFSS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	40.0 mm	1
60PFSS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	41.0 mm	1
60PFSS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	42.0 mm	1
60PFSS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	44.0 mm	1
60PFSS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	46.0 mm	1

5

Bivona® FlexTend™ Plus (Extra Length) Uncuffed Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFP40	4.0 mm	4.0 mm	6.0 mm	30.0 mm	44.0 mm	1
60PFP45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	48.0 mm	1
60PFP50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	50.0 mm	1
60PFP55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	52.0 mm	1
60PFP60	6.0 mm	6.0 mm	8.7 mm	30.0 mm	56.0 mm	1

6

Bivona® FlexTend™ Plus (Extra Length) Paediatric Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFPS40	4.0 mm	4.0 mm	6.0 mm	30.0 mm	44.0 mm	1
60PFPS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	48.0 mm	1
60PFPS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	50.0 mm	1
60PFPS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	52.0 mm	1
60PFPS60	6.0 mm	6.0 mm	8.7 mm	30.0 mm	56.0 mm	1

Bivona®

FlexTend™ TTS™

Reinforced proximal and distal sections provide a high level of kink and crush resistance. Extended tubes may reduce the risk of tracheostomy pressure ulcers on “no-neck” patients and deliver the advantages of a cuffed tube that acts like a non-cuffed tube when required.

1



67NFP30

2



67NFPS30

3



67PFS45

4



67PFSS45

5





67PFS45

6



67PFPS45



= Recommended for hospital use



= Recommended for home use



1

Bivona FlexTend™ with TTS™ Cuff Neonatal “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67NFP25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	30.0 mm	1
67NFP30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	32.0 mm	1
67NFP35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	34.0 mm	1
67NFP40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	36.0 mm	1

2

Bivona FlexTend™ with TTS™ Cuff Neonatal Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67NFPS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	30.0 mm	1
67NFPS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	32.0 mm	1
67NFPS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	34.0 mm	1
67NFPS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	36.0 mm	1

3

Bivona FlexTend™ with TTS™ Cuff Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67PFS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	38.0 mm	1
67PFS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	39.0 mm	1
67PFS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	40.0 mm	1
67PFS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	41.0 mm	1
67PFS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	42.0 mm	1
67PFS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	44.0 mm	1
67PFS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	46.0 mm	1

4

Bivona FlexTend™ with TTS™ Cuff Paediatric Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67PFSS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	38.0 mm	1
67PFSS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	39.0 mm	1
67PFSS35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	40.0 mm	1
67PFSS40	4.0 mm	4.0 mm	6.0 mm	20.0 mm	41.0 mm	1
67PFSS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	42.0 mm	1
67PFSS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	44.0 mm	1
67PFSS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	46.0 mm	1

5

Bivona FlexTend™ Plus (Extra Length) with TTS™ Cuff Paediatric “V” Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67PPF40	4.0 mm	4.0 mm	6.0 mm	30.0 mm	44.0 mm	1
67PPF45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	48.0 mm	1
67PPF50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	50.0 mm	1
67PPF55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	52.0 mm	1
67PPF60	6.0 mm	6.0 mm	8.7 mm	30.0 mm	56.0 mm	1

6

Bivona FlexTend™ Plus (Extra Length) with TTS™ Cuff Paediatric Straight Neck Flange Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
67PFPS40	4.0 mm	4.0 mm	6.0 mm	30.0 mm	44.0 mm	1
67PFPS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	48.0 mm	1
67PFPS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	50.0 mm	1
67PFPS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	52.0 mm	1
67PFPS60	6.0 mm	6.0 mm	8.7 mm	30.0 mm	56.0 mm	1




Bivona®

Sleep Apnea & Laryngectomy Tubes

The Bivona Sleep Apnea Tracheostomy Tube provides a low-profile design and is easily 'capped' and concealed during waking hours. The tube is intended for use in providing a patent airway in patients with a confirmed history of obstructive sleep apnea who do not require ventilator support.

The laryngectomy tube is designed in soft and compliant silicone. Available with wider neck flanges and longer shaft lengths.

The flexible all-silicone construction helps make both tubes comfortable.



= Recommended for hospital use



= Recommended for home use



1

Bivona® Sleep Apnea Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D	TUBE LENGTH	UNITS PER CASE
SAT605	6.0 mm	6.0 mm	9.0 mm	68 mm	1
SAT656	6.5 mm	6.5 mm	10.0 mm	73 mm	1
SAT757	7.5mm	7.5 mm	11.0 mm	79 mm	1

2

Bivona® Laryngectomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D	TUBE LENGTH	UNITS PER CASE
BOSL1L	9.5 mm	9.5 mm	12.0 mm	85 mm	1
BOSL1S	9.5 mm	9.5 mm	12.0 mm	55 mm	1
BOSL3L	11.0 mm	11.0 mm	14.0 mm	105 mm	1
BOSL3S	11.0 mm	11.0 mm	14.0 mm	55 mm	1
BOSL5L	13.0 mm	13.0 mm	16.0 mm	85 mm	1
BOSL5S	13.0 mm	13.0mm	16.0 mm	55 mm	1

Bivona[®]

Epistaxis Catheter

The Bivona[®] Epistaxis catheter offers rapid and positive control of nasal bleeding. The catheter design with its independently inflatable cuffs facilitates proper placement. The all silicone construction ensures that the cuffs remain soft and compliant, easily conforming to variations in nasal anatomy.

1



EP0050

PRODUCT CODE	SIZE I.D	FR	LENGTH	UNITS PER CASE
EP0050	5.0 mm	22	12.0 cm	1



= Recommended for hospital use

Bivona®

Customised

The extensive Portex® Bivona® Tracheostomy Tube range has a new offering: the opportunity to combine the benefits of the Bivona® silicone tube within your own customised design.

We recognise that you want to choose the most effective and comfortable tracheostomy option for every one of your patients, and this can sometimes mean you need to take the creative approach.

We offer a comprehensive customised tracheostomy tube service, so you can build a tube to meet your patients needs. A customised tube, designed by you, is a completely unique solution.

Our broad range of customisation options:

Tube Shaft Style

- Standard Silicone
- Hyperflex™ Wire Reinforced Silicone

Curvature

- Standard Anatomical
- Pre-curve Shaft
- Hyperflex™ Wire Reinforced Flexible

Length

- Standard Extensive Range

I.D. / O.D.

- Comprehensive range to choose from

Cuff Design


- Fome-Cuf® Cuff
- Aire-Cuf® Cuff
- TTS™ Cuff
- Cuffless

Neck Flange

- Adult
- Neonatal ("V" and "Straight")
- Paediatric ("V" and "Straight")
- Fixed
- Adjustable

Delivery options:

- Standard Service (non-sterile), 7 days
- Express Service (non-sterile), 3 days
- Same Day Service (non-sterile), 1-2 days (Only available in the US)
- Sterile Service, 4 weeks



Accessories

Code or Codes	Description	Unit
BRCA70	7.0 mm Bivona® Replacement Adult Inner Cannula	10
BRCA75	7.5mm Bivona® Replacement Adult Inner Cannula	10
BRCA80	8.0mm Bivona® Replacement Adult Inner Cannula	10
BRCA85	8.5 mm Bivona® Replacement Adult Inner Cannula	10
BRCA90	9.0 mm Bivona® Replacement Adult Inner Cannula	10
BRC270	7.0 mm Bivona® Replacement Adult Inner Cannula	2
BRC275	7.5 mm Bivona® Replacement Adult Inner Cannula	2
BRC280	8.0 mm Bivona® Replacement Adult Inner Cannula	2
BRC285	8.5 mm Bivona® Replacement Adult Inner Cannula	2
BRC290	9.0 mm Bivona® Replacement Adult Inner Cannula	2
BICSWB	Bivona® Inner Cannula Cleaning Swab	30
BV0818	Size 8 Tracheostoma Vent, 18.0 mm long	1
BV0827	Size 8 Tracheostoma Vent, 27.0 mm long	1
BV0836	Size 8 Tracheostoma Vent, 36.0 mm long	1
BV0855	Size 8 Tracheostoma Vent, 55.0 mm long	1
BV0918	Size 9 Tracheostoma Vent, 18.0 mm long	1
BV0927	Size 9 Tracheostoma Vent, 27.0 mm long	1
BV0936	Size 9 Tracheostoma Vent, 36.0 mm long	1
BV0955	Size 9 Tracheostoma Vent, 55.0 mm long	1
BV1018	Size 10 Tracheostoma Vent, 18.0 mm long	1
BV1027	Size 10 Tracheostoma Vent, 27.0 mm long	1
BV1036	Size 10 Tracheostoma Vent, 36.0 mm long	1
BV1055	Size 10 Tracheostoma Vent, 55.0 mm long	1
BV1227	Size 12 Tracheostoma Vent, 27.0 mm long	1
BV1236	Size 12 Tracheostoma Vent, 36.0 mm long	1
BV1255	Size 12 Tracheostoma Vent, 55.0 mm long	1
REDCAP	15mm Bivona® Decannulation Cap	26
CARSWL	Carden Swivel Connector	1
CON155	SidePort™ AutoControl™ Inflation Line Connector	50
CMD006	Cuff Maintenance Device With 60 cc Syringe & Stopcock	6
DRYPORT	DryPort Kit™ Inflation Line Connection Kit	1
DWG300	15mm Disconnection Wedge	25
TWILL	Twill Trach Tie	10



BRC290



DRYPORT



REDCAP



CMD006

Index

Symbols

60A150	15	60SP055	15	67PFSS55	23	BOSL5S	25
60A160	15	65N025	11	67SN025	9	BRC270	28
60A170	15	65N030	11	67SN030	9	BRC275	28
60A180	15	65N035	11	67SN035	9	BRC280	28
60A190	15	65N040	11	67SN040	9	BRC285	28
60AFHXL60	19	65P025	11	67SP025	9	BRC290	28
60AFHXL70	19	65P030	11	67SP030	9	BRCA70	28
60AFHXL80	19	65P035	11	67SP035	9	BRCA75	28
60AFHXL90	19	65P040	11	67SP040	9	BRCA80	28
60HA25	17	65P045	11	67SP045	9	BRCA85	28
60HA30	17	65P050	11	67SP050	9	BRCA90	28
60HA35	17	65P055	11	67SP055	9		
60HA40	17	65SN025	11	75FHXL60	19	C	
60HA45	17	65SN030	11	75FHXL70	19	CARSWL	28
60HA50	17	65SN035	11	75FHXL80	19	CMD006	28
60HA55	17	65SN040	11	75FHXL90	19	CON15S	28
60NFP25	21	65SP025	11	75HA60	17		
60NFP30	21	65SP030	11	75HA70	17	D	
60NFP35	21	65SP035	11	75HA80	17	DRYPORT	28
60NFP40	21	65SP040	11	75HA90	17	DWG300	28
60NFPS25	21	65SP045	11	85N030	13		
60NFPS30	21	65SP050	11	85N035	13	E	
60NFPS35	21	65SP055	11	85N040	13	EP0050	26
60NFPS40	21	67FHXL60	19	85P025	13		
60N025	15	67FHXL70	19	85P030	13	R	
60N030	15	67FHXL80	19	85P035	13	REDCAP	28
60N035	15	67FHXL90	19	85P040	13		
60N040	15	67HA60	17	85P045	13	S	
60P025	15	67HA70	17	85P050	13	SAT605	25
60P030	15	67HA80	17	85P055	13	SAT656	25
60P035	15	67HA90	17	670150	9	SAT757	25
60P040	15	67N030	9	670160	9		
60P045	15	67N035	9	670170	9	T	
60P050	15	67N040	9	670175	9	TWILL	28
60P055	15	67NFP25	23	670180	9		
60PFP40	21	67NFP30	23	670185	9		
60PFP45	21	67NFP35	23	670190	9		
60PFP50	21	67NFP40	23	750150	11		
60PFP55	21	67N025	9	750160	11		
60PFP60	21	67P025	9	750170	11		
60PFPS40	21	67P030	9	750180	11		
60PFPS45	21	67P035	9	750190	11		
60PFPS50	21	67P040	9	750195	11		
60PFPS55	21	67P045	9	755150	11		
60PFPS60	21	67P050	9	755160	11		
60PFS25	21	67P055	9	755170	11		
60PFS30	21	67PFP40	23	755180	11		
60PFS35	21	67PFP45	23	755190	11		
60PFS40	21	67PFP50	23	755195	11		
60PFS45	21	67PFP55	23	850150	13		
60PFS50	21	67PFP60	23	850160	13		
60PFS55	21	67PFPS40	23	850170	13		
60PFSS25	21	67PFPS45	23	850180	13		
60PFSS30	21	67PFPS50	23	850190	13		
60PFSS35	21	67PFPS55	23	850195	13		
60PFSS40	21	67PFPS60	23	855150	13		
60PFSS45	21	67PFS25	23	855160	13		
60PFSS50	21	67PFS30	23	855170	13		
60PFSS55	21	67PFS35	23	855180	13		
60SN025	15	67PFS40	23	855190	13		
60SN030	15	67PFS45	23			B	
60SN035	15	67PFS50	23	BICSWB	28	BOSL1L	25
60SN040	15	67PFS55	23	BOSL1S	25	BOSL3L	25
60SP025	15	67PFSS25	23	BOSL3S	25	BOSL3S	25
60SP030	15	67PFSS30	23	BOSL5L	25		
60SP035	15	67PFSS35	23				
60SP040	15	67PFSS40	23				
60SP045	15	67PFSS45	23				
60SP050	15	67PFSS50	23				



PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical Sales Office
Country Sales Office Address
Phone: Sales Office
Country email or Web Address

Find your local
contact information at:
[www.smiths-medical.com/
customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Portex, Bivona, TTS, Aire-Cuf, DryPort Kit, Fome-Cuf, Hyperflex, FlexTend and the Portex and Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2014 Smiths Medical. All rights reserved. MHYTCA-1053. TR194451GL-012014

smiths medical